# LUBEXXX®

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**UK Instructions for use**

**LUBEXXX® Pelvic Floor Trainer** for stronger pelvic floor muscles is available in two different weights: 72g and 92g for a moderately trained pelvic floor.

Both trainer stimulate, massage and train your muscles. The gentle to intense vibrations of the trainer give your pelvic floor muscles a perfect workout – whenever and wherever you want. A well-trained pelvic floor enhances a woman's ability to reach orgasm, making sex more pleasurable. The trainer is used vaginally like a conventional tampon. When inserted, the tear-resistant retrieval string should be easy to reach. Please relax your pelvic floor during removal. For easier insertion and removal, please use LUBEXXX® Original lubricant.

**Use:** To strengthen the pelvic floor muscles the pelvic floor trainer is inserted lengthwise into the vagina, similarly to a tampon, with the retrieval string remaining outside. The vibrations of the inserted balls send a signal to the pelvic floor that there is something it needs to hold on to, thereby strengthening its muscles. In average a daily training of 8 to 10 minutes is recommended. For best results, they can be worn during everyday activities. They are pulled back out using the tear-resistant retrieval string. For more information on pelvic floor training, please ask your physician. Women who would like to train their pelvic floor during pregnancy should only do this after consulting their midwife or physician. After giving birth, pelvic floor trainer should only be used after any birth-related injuries have healed, the lochia has stopped and postnatal exercises have been started.

**Care:** Before first use and after each use, the pelvic floor trainer should be cleaned with LUBEXXX® toy cleaner or using lukewarm water and mild soap.

**Warning:** The pelvic floor trainer is not suitable for anal use. Please remove the trainer before vaginal intercourse. The two pelvic floor exercisers are not a children's toy and should be kept away from children.

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MAKE HoBo marketing GmbH,